



Remote health monitoring



PhD Candidate Sonja Babac

13th Jan 2025

About me



Present:

- PhD Candidate @ TU/e – Philips – Catharina Hospital
- Researching AI & Signal Processing for Remote Health Monitoring


Past:

- Graduated Electrical Engineering Master @ TU/e in 2023
- Graduate Intern @ Philips - Vital Signs Camera
- Intern @ Nemo Healthcare - Remote Fetal Monitoring

What is meant by remote health monitoring?



Health Continuum



Remote Health Monitoring

Clinical Home Care:


- Patient monitoring using **clinical-grade** wearables
- Continuous postoperative vital signs data
- Telemonitoring for chronic diseases
- Telemonitoring for outpatient clinics
- Telehealth consultations

Fitness & prevention:

- health monitoring using **fitness** wearables
 - Cardiovascular health
 - Sleep
 - Recovery
 - Stress
 - Activity
 - Women's Health
 - Metabolic Health
- Health data summaries for physicians

What are the advantages
of remote patient
monitoring?





Remote Health Monitoring in Clinical Context

Remote patient monitoring

“The use of a device for interaction between providers and patients outside of the provider’s organization “



Why remote patient monitoring?

“

Many patients want to recover in their own environment. Supported by the Healthdot, we offer these people the option to go home earlier with the guarantee that the hospital is still monitoring their health. This form of continuous monitoring also offers the possibility to detect complications at an early stage. The hospital can therefore provide good care while relieving the burden on bed capacity.”

Dr. Simon Nienhuijs
Bariatric Surgeon at Catharina Hospital



73%

of healthcare professionals across 19 countries say connected technology is important to improving the prevention of medical issues.²²



By 2022,

60%

of chronic disease pathways will involve remote patient monitoring.¹⁶



30%

of all ambulatory encounters will be virtual, up from about 2% at the start of 2020.¹⁶

Why remote patient monitoring?

- Data-driven clinical decision making
 - Accurate vitals signs monitoring guaranteed
 - High-frequency continuous measurements
 - Early deterioration detection
- Lowering healthcare costs
 - Relieving the burden on bed capacity
 - Relieving the burden of staff shortages
- Home comfort
 - Receive health guidance @ home remotely (Telehealth)
 - Fewer commutes to the hospital
 - Lowering chances of hospital infections
 - Recover or monitor health in your own environment

Postoperative remote monitoring



Postoperative remote monitoring



Patient receives Healthdot sensor



Care provider can monitor patient at home



If there is anything, the healthcare provider can respond


Early deterioration detection

- Cardiovascular health
 - Heart rate
 - Resting heart rate
 - Heart rate variability
 - Respiration rate
 - Cardiorespiratory coupling
- Sleep health
 - Circadian rhythm
 - Sleep quality
- Activity
 - Daily movement



Early warning score

Remote Early Warning Score




The REWS was built based solely on automatically collected HR and RR measurements from the wearable.

Receiver operating characteristic (ROC) curve for the Modified Early Warning Score (MEWS) and Remote Early Warning Score (REWS) to describe the development of a deterioration in the next 24 h. Displaced numbers are the area under the ROC-curves. Dashed line represents the random-guessing line.

[van der Stam, J. A. et al.](#)


Circadian rhythm in postoperative patients



Postoperative heart rate of the bariatric (blue) and cancer surgery (green) subgroups. The line indicates the median heart rate, and the shaded area indicates the IQR.

van der Stam, J. A. et al.

Circadian rhythm in postoperative patients



Postoperative peak, nadir and peak-nadir excursions of the heart rate of the bariatric (blue) and cancer surgery (green) subgroups. Dots indicate the daily median value for the respective features, and shaded areas indicate the IQR.

van der Stam, J. A. et al.

Devices for postoperative patient monitoring



Health Dot

Accelerometry based:

- Heart rate
- Respiration rate
- Activity
- Posture



[Video link](#)

ViQtor

PPG based:

- Heart Rate
- HRV
- Respiration Rate
- SpO2

Accelerometry based:

- Activity

Other:


- Temperature



Other use cases of remote patient monitoring



Clinical grade sleep monitoring @ home




**Sleep diagnostics
and monitoring,
redefined.**

The first self-applied, no-wire
Polysomnography solution.




[Video link](#)

Clinical grade fetal monitoring @ home



Tracing with the Nemo Fetal Monitoring System



Standard monitoring with fetal scalp electrode (FSE) and toco

"We need a non-invasive fetal monitoring system - also in a future with obesity - which provides qualitative information and allows intervention, before any damage occurs to the unborn child."

Prof. Dr. Y. (Yves) Jacquemyn, UZA - Antwerp University Hospital.



[Video link](#)


Chronic disease monitoring @ home

- COPD
- Asthma
- Pacemaker patients
- Heart Failure
- Atrial Fibrillation
- Diabetes
- Obesity
- Hypertension



[Video link](#)

Chronic disease monitoring @ home



Chronic disease monitoring @ home



Diabetes monitoring @ home


- Continuous Glucose Monitor
- App for diabetes management
 - Food intake
 - Insulin intake
 - Medication management
 - Alerts for low blood glucose
 - Alerts for high blood glucose



Vital Signs Camera



Vital Signs Camera




[Video link](#)

Vital Signs Camera

- Uses remote PPG (PhotoPlethysmoGraphy)
- Camera tracks rise/fall of chest and abdomen to determine respiration rate
- Tiny changes in skin color caused by HR are detected (invisible micro-blushes)
- Philips solved many challenges associated with this technique (motion, light, etc.)
- Suitable for TeleHealth consultations
- Suitable for monitoring @ lower acuity wards

Vital Signs Camera




[Video link](#)

Vital Signs Camera



[Video link](#)

Vital Signs Camera



[Video link](#)

Remote Health Monitoring for Fitness & Prevention





Fitness wearables: wrist band




[Video link](#)

Fitness wearables: smart ring



24/7 Heart Rate Monitoring

See how your heart rate changes, day and night. Learn how your body responds to daily habits, meals, stress, and environment.




Cardiovascular Age

Quickly understand how your cardiovascular system is aging relative to your chronological age, and what that may mean for long-term health.



Cardio Capacity

Get an estimate of your VO2 Max, or how efficiently your body delivers and utilizes oxygen during intense physical effort.



Heart Rate Variability

Oura tracks your HRV to help you understand if you're stressed, sick, or need to take time to recover.

Fitness wearables: smart ring




Sleep Score

The answer to "how did you sleep last night?" Your score is based on your total sleep, heart rate variability (HRV), nighttime movement, sleep regularity, and more.




Sleep Stages

Wake up to in-depth analysis of your deep sleep, REM sleep, and light sleep from the night prior.



Blood Oxygen Sensing


By detecting your blood oxygen levels at night, Oura can tell if you're experiencing any breathing disturbances.



Bedtime Guidance

Oura reminds you when to start winding down so you can improve your sleep quality.

Women's Health Monitoring



Health monitoring in-vehicles



[Video link](#)

Adoption outlook of remote health monitoring





CONNECTED AT HOME



Interested in an MSc internship or
thesis project in Remote Health
Monitoring?

Contact: s.babac@tue.nl



Thank you for your attention

e/MTIC contact:
info@emtic.nl | www.emtic.nl  e-mtic

e/MTIC Partners:



Kempenhaghe



máxima
medisch centrum

PHILIPS

TU/e EINDHOVEN
UNIVERSITY OF
TECHNOLOGY